



Application

SECUTAPE® fixing system GECKO

GECKO holds reliably, without any adhesive tubes and wires with the new patented **NANOPLAST®** antislip foil.

The fixation can be **opened and closed many times** without leaving any marks on the tube!





Gecko comes as a set with different holders, also in Hydrocolloid.

The fixation **stays even when it gets wet.**

NANOPLAST® has the so-called Gecko effect and is made out of silicon with aprox. 29000 nano legs per qcm.

Pull the Velcro fastener through, **the green side** towards the hose. Tighten slightly with both hands and place on the lower part.

	Velcro-Binder ord.no. 8501000 PU 200 St.
	Velcro-binder - Set SMAL Hydrocolloid ord.no. 8501100 PU 100 St.
	Velcro-binder - Set BIG ord.no. 8501130 PU 100 St.



TechniMed AG

SWITZERLAND since 1987

fon 0041 / 71 844 49 20

Kirchstrasse 54

CH-9400 Rorschach

fax 0041 / 71 845 52 45

Germany

fon 0800 007 9892

fax 0800 941 3012

Austria

fon 01/ 440 7979

fax 01/ 440 7990



Made in Germany

www.secutape.ch